

Champion E-33-12V

Fixed Angle Portable Centrifuge

Centrifuge

www.AmpleScientific.com

(678)380-0369



Features

- Portable design with handle for secure transportation
- Small footprint with suction feet for stability
- Standard 12V import voltage to be used on automobiles
- Optional ACDC adaptor allows centrifuge to be used indoors
- Auto-stop system – automatically decelerate when lid is unlocked
- Additional set (8) of small tube sleeves to accommodate 7ml tubes

Specifications

Item	Specification
Speed	3,300 rpm (fixed Speed)
Max Rcf	1350 xg
Timer	Analog, 0 – 30 min
Noise at Max Speed	≤65dB
Power	12V automotive power
Motor	Brushless
Weight	10 lb
Dimensions (LxWxH)	320 x 320 x 320 mm
Rotor	45° angle 8 place PC rotor With 8x15ml PC buckets
Compatible Test Tubes	15ml test tubes w/ diameter less than 18mm 7ml test tubes w/diameter less than 13mm