

Champion E-33

Fixed Angle Bench-top Centrifuge

Centrifuge

www.AmpleScientific.com



Features

- Small footprint with suction feet for stability
- Auto-stop system – automatically decelerates when lid is unlocked
- Temperature controlled brushless motor
- Spins up to 8 of 15ml test tubes
- Additional set (8) of small tube sleeves to accommodate 7ml tubes

Specifications

Item	Specification
Speed	3,300 rpm (fixed Speed)
Max Rcf	1350 xg
Timer	Analog, 0 – 30 min
Noise at Max Speed	≤65dB
Power	110V 60Hz or 220V 50Hz
Motor	Brushless
Weight	10 lb
Dimensions (LxWxH)	320 x 320 x 320 mm
Rotor	45° angle 8 place PC rotor With 8x15ml PC buckets
Compatible Test Tubes	15ml test tubes w/ diameter less than 18mm 7ml test tubes w/diameter less than 13mm