

# Champion F-33D

Fixed Angle Bench-top Centrifuge

Centrifuge

www.AmpleScientific.com

(678)380-0369



## Features

- Vortex ventilation system keeps sample temperature consistent, while maintaining sample integrity
- Wide mouth opening allows easy access to inner chamber
- Closed system inner chamber facilitates cleaning and maintenance
- Automatically stops acceleration when lid opens
- Temperature controlled maintenance -free brushless motor
- Digital tachometer and RPM adjustment allows precise speed control
- Additional set (8) of small tube sleeves to accommodate 7ml tubes

## Specifications

Item	Specification
Speed	500 - 3,300 rpm
Max Rcf	1350 xg
Timer	Digital, 0 – 30 min
Noise at Max Speed	≤65dB
Power	120V 60Hz or 220V 50Hz
Motor	Brushless
Weight	18 lb
Dimensions (LxWxH)	400 x 320 x 250 mm
Rotor	45° angle 8 place PC rotor With 8x15ml PC buckets
Compatible Test Tubes	15ml test tubes w/ diameter less than 18mm 7ml test tubes w/diameter less than 13mm