

Champion S-33

Swing Out Rotor Bench-top Centrifuge

Centrifuge

www.AmpleScientific.com

(678)380-0369



Standard S33-SR8 Rotor



S33-SR4 Rotor

Features

- Vortex ventilation system keeps sample temperature consistent, while maintaining sample integrity
- Wide mouth opening allows easy access to inner chamber
- Closed system inner chamber facilitates cleaning and maintenance
- Temperature controlled maintenance-free brushless motor
- Digital tachometer and RPM adjustment allows precise speed control
- Smooth acceleration and deceleration
- Optional rotor for 50ml tubes

Specifications

| Item | Specification |
|-----------------------|--|
| Speed | 300 - 3,300 rpm |
| Max Rcf | 1732 xg |
| Timer | Digital, 0 – 30 min |
| Noise at Max Speed | ≤65dB |
| Power | 110V 60Hz or 220V 50Hz |
| Motor | Brushless |
| Weight | 30 lb |
| Dimensions (LxWxH) | 462 x 362 x 318 mm |
| Rotor | Standard [S33-SR8] : 8x15ml duralumin swing out rotor with PC buckets <i>Optional [S33-SR4] :</i> 4x50ml duralumin swing out rotor with PC buckets |
| Compatible Test Tubes | 15ml test tubes w/ diameter <18mm 50ml test tubes w/ diameter <30mm |